
A

Card A: Preparing Culture Media for 10 Tests.

A1

Boil water for 10 minutes, allow to cool.

A2

Add 10 spoonfuls of culture media into measuring tube using the integral spatula.

A3

Fill the Media Measuring Device with the cooled sterilised water to the bottom lip.

A4


Close the lid, and then shake to dissolve all the powder.

A5

The dissolved culture medium should turn a bright pink colour.

A6

Sterilise the Media Measuring Device at 121°C for at least 10 minutes.

B

Card B: Sterilising the Membrane Filtration Unit.

B1

Add 1ml of methanol into the sampling cup, swirling to coat internal surfaces.

B2

Push the filter funnel into the blue rubber base leaving sterilisation symbol visible.

B3

Hold the sampling cup to face away from you. Use a lighter to ignite the methanol.

B4

When flame is nearly extinguished, invert funnel and rubber insert into the cup.

B5

Wait 15 minutes.

B6

Rinse the sterile sampling cup with sample water before use.

C

Card C: Preparation of the Aluminium Petri Dishes and Absorbent Pads.

C1

Remove the lid from the tube of absorbent pads and attach to the pad dispenser.

C2

Remove lid of sterile petri dish taking care only to handle the lid by the sides.

C3

Dispense absorbent pad into petri dish by sliding the grooved lever of the dispenser.

C4

Carefully pour the media from the MMD on to the absorbent pad in the petri dish.

C5

Ensure the absorbent pad is saturated with a small excess of media at the edges.

C6

Replace petri dish lid to keep contents sterile and free from contamination.

D

Card D: Membrane Filtration of the Sample.

D1

Pull filter funnel from rubber base support. Invert filter funnel, place on clean work-surface.

D2

Sterilise forceps tips by passing them through flame from a lighter. Allow them to cool.

D3

Peel back clear membrane filter outer wrapper. Use sterile forceps to separate the white, gridded filter from blue backing paper.

D4

Place the membrane filter directly onto the bronze filter support disc ensuring the gridded side is face-up.

D5

Lock membrane filter in place by pushing filter funnel firmly down into position in blue rubber base. Ensure filter funnel is in the "Filtration" position.

D6

Pour water sample into filter funnel up to 100ml graduation. Connect hand-vacuum pump to membrane filtration unit and use pump to filter the sample.

D7

When water is filtered, remove filter funnel & vacuum pump from rubber base. Remove membrane filter with sterile forceps; place it grid-side-up on wetted growth pad in petri dish.

D8

Replace the petri dish lid and label appropriately. Place into the petri dish rack until you are ready to incubate.

E

Card E: Counting Colonies and Recording Results.

E1

After incubation remove petri dishes.
Note incubation temperature and period.

E2


Place petri dishes on a flat surface in
good natural light. Avoid direct sunlight.

E3


Remove lids. Count yellow colonies more
than 1mm diameter. Use lens if necessary.

E4


DO NOT count colonies that are
CLEAR, RED or ANY OTHER COLOUR.

E5

Use the grid on the paper to help count
colonies systematically.

E6

= 13 CFU/100ml.

Results expressed as colony forming
units per unit volume; CFU/100ml.

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